Orbegozo

ALMOHADILLA ELÉCTRICA – MANUAL DE INSTRUCCIONES
HEATING PAD - INSTRUCTION MANUAL
COUSSIN CHAUFFANT - MANUEL D'INSTRUCTIONS
ALMOFADA ELECTRICA - MANUAL DE INSTRUÇÕES

AH 3050 (40X30CM) AH 3150 (50X38CM) AH 3160 (76X43CM)

Sonifer, S.A.
Avenida de Santiago, 86
30007 Murcia
España
E-mail: sonifer@sonifer.es
Made in P.R.C.

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GENERAL SAFETY INSTRUCTIONS

1. This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children must never play with the appliance. Cleaning and user maintenance must never be carried out by children without supervision.
2. Children should be supervised to ensure that they do not play with the appliance.
3. Keep the appliance and the cable cord out of the reach of children less than 8 years old.
4. **WARNING:** In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach.
5. If the supply cord is damaged it must be repaired by the Authorized Service Agent, in order to avoid hazards.
6. Never pull on the cord when unplugging.
7. Do not use the unit with a damaged cord or plug, or if it is not working properly.
8. Do not handle the appliance with wet hands.
9. Never immerse the appliance in water or any other liquid.
10. Make sure the appliance has been unplugged before cleaning.
11. This appliance must be installed following the national regulations for electrical installations.
12. This appliance is for household use only.

13. In case that you need a copy of the instruction manual, you can ask for it writing an email to sonifer@sonifer.es.
14. WARNING: In case of misuse, there is a risk of possible injury.
- Important instructions. Retain for future.
 - During use the user must not fall asleep.
 - When the device is not operating, disconnect it. The heating blanket should be stored in a dry place.
 - Do not use the product when it is folded, and do not iron out creases.
 - Inspect the product regularly for signs of wear or damage. If there are signs of wear or damage or if the product has been used incorrectly, return it to the distributor before further use.
 - This product is not designed for medical use in hospitals.
 - The pad should be used under direct supervision, especially when used by people who are disabled, insensitive to heat or sick or when used by children.
 - The cover can be washed. The rest should not be immersed in water or other liquid.
 - If the device is used for a prolonged period on the highest setting, burning of the skin may result.
 - If the product is used by children, the thermostat should be set by a parent or supervisor. Use of the device must be monitored!
 - Do not use the product if it is wet or in rooms with a high level of humidity, such as bathrooms etc.
 - The heating pad must always be used inside the cover.

- To put it into operation, put the regulator in position for maximum 20 minutes to obtain a rapid heating and then switch to the low position. The position of maximum heat is used only for the rapid warming and never more than 20 minutes.
- Do not fold or bend the electric cord and do not use your blanket when it is folded up, rolled up or bent on itself.
- Wear it over a garment.
- Do not cut or pierce the heating blanket. Do not make holes using pins or dull instruments. Do not disassemble it
- Do not wash or dry clean the inner part of the heating blanket.

OPERATION

The heating pad can smell plastic during the first use but it will disappear after a short period of time.

ON/OFF

To use the device, first connect the controller to the pad and then plug the controller into a power outlet.

Switch on:

Select the highest temperature level (level 6) a short period of time before use, so the pad will warm up faster. Select the level of temperature control 1 to 6. If you use the pad for a few hours select an average temperature (level 3).

Shutdown:

Select the level 0 on the remote to turn off the heating pad and unplug from the wall outlet.

TEMPERATURE LEVELS

Level 0: off.

Level 1-2: minimum heat.

Level 3-4: medium heat.

Level 5-6: maximum heat.

AUTO STOP

The heating pad is turned off automatically after 90 minutes. It will appear "A" in the LED indicator on the controller.

To reuse the pad, first select the 0 position on the knob and after about 5 seconds select the desired level to restart the pad.

If you do not want to continue using the pad after auto stop, select level 0 in the controller and unplug it from the outlet.

CLEANING AND CARE

Use the heating pad always with the fabric cover to prevent soiling or damage to the appliance (particularly if you use essential oils or creams).

Cleaning: the pad cover is washable. The washing method is specific.

Examine occasionally the pad to verify that the flexible cord shows no traces of wear and the regulator is not damaged. If you find any damage discontinue use of the pad and take it for service.

STORAGE

Wait for the heating pad to cool before folding and storing, otherwise it could result in damaged to the appliance.

To prevent the pad is damaged, do not leave any objects on it when storing.

If you are not using the heating pad for a long time keep it in its original packaging in a dry place and no heavy objects on top.

MEANING OF SYMBOLS APPEARING ON HEATING PAD

Please, read de instructions.

Do not use if folded is creased.

Do not insert pins or needles.

Do not use bleach.

Wash up 30°C

The cover can be ironed

Do not dry clean

Do not dry the cover or blanket in the dryer.

Do not use on babies.

The blanket cannot be ironed.

Blanket cannot be washed.

Disposal of old electrical appliances.

The European directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE), requires that old household electrical appliances must not be disposed of in the normal unsorted municipal waste stream. Old appliances must be collected separately in order to optimize the recovery and recycling of the materials they contain, and reduce the impact on human health and the environment.

The crossed out "wheeled bin" symbol on the product reminds you of your obligation, that when you dispose of the appliance, it must be separately collected.

Consumers should contact their local authority or retailer for information concerning the correct disposal of their old appliance.

DECLARATION OF CONFORMITY:

This device complies with the requirements of the Low Voltage Directive 2006/95/EC and the requirements of the EMC directive 2004/108/EC.

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